



# Hangzhou New Asia International Co., Ltd

Tel:86-571-85815391 Fax:86-571-85815449

Email:Sales@herbalsextract.com

---

## Determination of $\beta$ 1-6 linked $\beta$ 1-3 glucans in the yeast cell walls (Wolff Method)

### 1 . Principle

The YCW is first submitted to hydrolysis with NaOH. The insoluble part is hydrolyzed by H<sub>2</sub>SO<sub>4</sub> and, after neutralization, the glucose content is determined by an enzymatic kit (Roche)

### 2 . Material

- 50 ml capped flask
- lab centrifuge
- centrifuge tube
- pH - meter
- filter GF/C Whatmann
- 500 ml flask
- oven
- NaOH 0,1 N
- NaOH 32 %
- H<sub>2</sub>SO<sub>4</sub> 72 %
- Enzymatic kit for glucose determination (Boehringer/Roche) ref : 716251

### 3 . Method

In a 50 ml capped flask

- ◆ Weigh 2 g of YCW in a 50 ml capped flask.
- ◆ note the exact weight W<sub>5</sub> and calculate the ratio  $RW_5 = 1/W_5$  for result in § 4
  
- ◆ Add 20 ml of 0,1 N, NaOH.
- ◆ Stir with magnetic stirring overnight at room temperature
- ◆ Put the solution into a centrifuge tube.
- ◆ Centrifuge 10 min
- ◆ Remove the supernatant and disperse the pellet in 40 ml of water
- ◆ Repeat 3 time this centrifuge operation
- ◆ Remove the supernatant after the last centrifuge
- ◆ In the centrifuge tube containing the pellet, add 25 ml of H<sub>2</sub>SO<sub>4</sub> 36 N.
- ◆ → Don't add the 25 ml of acid in one time. Add the acid slowly (5 ml by 5 ml) and stir the tube after each addition.
- ◆ Stir during 30 min at room temperature
- ◆ Put the entire solution into a capped flask and add 300 ml of water
- ◆ Close the flask and leave at 100°C during 6 H in a oven
- ◆ After cooling, neutralize at pH 4,5 with NaOH.
- ◆ Complete the solution at 500 ml.
- ◆ Filter with Wathman GFA membrane
- ◆ Determine the glucose in the filtrate with enzymatic kit (Roche) (see § 3.3.)